

Portable Proteomics Pipeline (P3)

Use Case: PrideID

Daniel Kristiyanto (daniel.kristiyanto@pnnl.gov)

June 24, 2016

This is a use case/manual for a P3 pipeline with files retrieved from Pride database.

1. Make sure Docker engine is installed and running.

```
docker ps
```

2. Download a sample p3.config and MSGFDB_Mods.txt file and put it to a directory. e.g.
~/Desktop/PXD001468

```
cd ~/Desktop/PXD001468
```

in this case, p3.config must have the following fields:

```
REPO = FTP
```

```
FTP1 = (e.g: ftp://ftp.pride.ebi.ac.uk/pride/data/archive/2015/06/PXD001468)
```

```
RUN_MSGF = (e.g: YES)
```

```
METHOD = (e.g: SPECTRUM_COUNT)
```

If other PrideID is desired, change the PRIDEID option in the p3.config. Ensure that the corresponding PrideID have the required files (*.fasta, *.mzml, etc.).

3. Download / update P3 from Dockerhub, and run the image.

```
eval $(docker-machine env default)
```

```
docker pull kristiyanto/p3
```

```
docker run --rm -v ~/Desktop/PXD001468:/root/data kristiyanto/p3
```

if necessary, another image can be run simultaneously (in a different terminal windows).

```
docker run --rm -v ~/Desktop/PXD001468:/root/data kristiyanto/p3
```

4. Once processes completed, the directory will be populated with:

- Mass Spectrometry files downloaded from Pride repository (*.mzid, *.mzml)
- MSGF+ output files (*.mzid, *.canon, *.revCat.cnlcp, *.revCat.csarr, *.revCat.cseq)
- Output table (*.txt)
- Output file as R objects (*.rda)

If the process stop returning error: * Make sure all the required files are in the folder. E.g: Fasta files sometimes are not provided in the FTP/Pride repository and must be provided manually. * Remove all *.tmp files before re-run the containers.