## Portable Proteomics Pipeline (P3) Use Case: PrideID

Daniel Kristiyanto (daniel.kristiyanto@pnnl.gov)

June 23, 2016

## Pride ID

This is a use case/manual for a P3 pipeline with files retrieved from Pride database.

1. Make sure Docker engine is installed and running.

docker ps

2. Download a sample p3.config and MSGFDB\_Mods.txt file and put it to a directory. e.g.

```
cd ~/Desktop/PXD001468
ls
```

in this case, p3.config must have the following fields:

```
REPO = PRIDE
PRIDEID = (e.g: PXD001468)
RUN_MSGF = (e.g: YES)
METHOD = (e.g: SPECTRUM_COUNT)
```

If other PrideID is desired, change the PRIDEID option on the p3.config. Ensure that the corresponding PrideID have the required files (\*.fasta, \*.mzml, etc.).

3. Download / update p3 from Dockerhub, and run the image.

```
docker pull kristiyanto/p3
docker run --rm -v ~/Desktop:/root/data kristiyanto/p3
```

if necessary, another image can be run simultaneously (in a different terminal windows).

- 4. Once processes completed, the directory will be populated with:
- Mass Spectometry files downloaded from Pride repository (\*.mzid, \*.mzml)
- Peptide Sequence file (\*.fasta)
- MSGF+ output files (\*.mzid, \*.canno, \*.revCat.cnlcp, \*.revCat.csarr, \*.revCat.cseq)
- Output table (\*.txt)
- Output file as R objects (\*.rda)